
AprioCore[®] *plus*

Semi-automatic biopsy system

- (EN) Instruction for use
- (SV) Bruksanvisning
- (NL) Handleiding
- (FR) Mode d'emploi
- (ES) Instrucciones de uso
- (IT) Istruzioni per l'uso
- (DE) Gebrauchsanweisung
- (NO) Bruksanvisning

AprioMed

(EN) Table 1. Device information (SV): Tabell 1. Produktinformation (NL) Tabel 1. Productinformatie (FR) Tableau 1. Informations sur le dispositif (ES): Tabla 1. Información del dispositivo (IT) Tabella 1. Informazioni sul dispositivo (DE) Tabelle 1. Instrumenteninformationen (NO) Tabell 1. Produktinformasjon

(EN) Gauge x Length (SV) Gauge x Längd (NL) Gauge x Lengte (FR) Jauge x Longueur (ES) Calibre x longitud (IT) Calibro x lunghezza (DE) Durchmesser x Länge (NO) Gauge x lengde	(EN) Specimen notch size (SV) Storlek skåra för vävnadsprov (NL) Lengte biopt (FR) Encoche de l'échantillon à prélever (ES) Tamaño de la muesca de la muestra (IT) Dimensioni del campione (DE) Größe der Biopsiekammer (NO) Størrelse av hakk for vevsprøve	(EN) Cannula throw length (SV) Slaglängd kanyl (NL) Penetratiediepte (FR) Canule Profondeur de pénétration (ES) Penetración de la cánula (IT) Lunghezza della cannula (DE) Größe der Wurfänge (NO) Kanylens slaglengde	(EN) Dead space length (SV) Längd distalt dödrymme (NL) Lengte loze ruimte (FR) Longueur de l'espace mort (ES) Longitud del espacio muerto (IT) Lunghezza dello spazio morto (DE) Länge des Totraums (NO) Lengde distalt dødrum	(EN) REF AprioCore plus (SV) REF AprioCore Plus (NL) REF AprioCore plus (FR) REF AprioCore plus (ES) REF AprioCore plus (IT) CODICE AprioCore plus (DE) REF AprioCore Plus (NO) REF AprioCore plus	(EN) REF compatible coaxial needle (SV) REF kompatibel koaxialnål (NL) REF bijpassende coaxiaalnaald (FR) REF Aiguille coaxiale compatible (ES) REF Aguja coaxial compatible (IT) CODICE Ago coassiale compatibile (DE) REF kompatibel Koaxialnadel (NO) REF kompatibel koaxialnål
16G x 8,5 cm	10 mm	16 mm	2,7 mm	<input type="radio"/> ACP16-085-10	SG15-071/HG15-071
16G x 13,5 cm	10 mm	16 mm	2,7 mm	<input type="radio"/> ACP16-135-10	SG15-121/HG15-121
16G x 18,5 cm	10 mm	16 mm	2,7 mm	<input type="radio"/> ACP16-185-10	SG15-171/HG15-171
16G x 9,5 cm	20 mm	26 mm	2,7 mm	<input type="radio"/> ACP16-095-20	SG15-071/HG15-071
16G x 14,5 cm	20 mm	26 mm	2,7 mm	<input type="radio"/> ACP16-145-20	SG15-121/HG15-121
16G x 19,5 cm	20 mm	26 mm	2,7 mm	<input type="radio"/> ACP16-195-20	SG15-171/HG15-171
18G x 8,5 cm	10 mm	16 mm	2 mm	<input checked="" type="radio"/> ACP18-085-10	SG17-071/HG17-071
18G x 13,5 cm	10 mm	16 mm	2 mm	<input checked="" type="radio"/> ACP18-135-10	SG17-121/HG17-121
18G x 18,5 cm	10 mm	16 mm	2 mm	<input checked="" type="radio"/> ACP18-185-10	SG17-171/HG17-171
18G x 9,5 cm	20 mm	26 mm	2 mm	<input checked="" type="radio"/> ACP18-095-20	SG17-071/HG17-071
18G x 14,5 cm	20 mm	26 mm	2 mm	<input checked="" type="radio"/> ACP18-145-20	SG17-121/HG17-121
18G x 19,5 cm	20 mm	26 mm	2 mm	<input checked="" type="radio"/> ACP18-195-20	SG17-171/HG17-171

(EN)

- 1 Trigger button
- 2 Specimen release button
- 3 Cutting cannula with centimeter marking (marked in 1 cm increments) and echogenic distal tip
- 4 Specimen notch
- 5 Cannula throw length
- 6 Distal dead space

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





ENGLISH

DEVICE DESCRIPTION

AprioCore plus is a sterile, disposable, semi-automatic, spring-loaded biopsy needle. AprioCore plus is provided in a variety of sizes and can be used together with coaxial needles such as AprioMed disposable Coaxial Needles Gangi-SoftGuard and Gangi-HydroGuard, see table 1.

INTENDED USE

AprioCore plus is intended for use in obtaining core biopsy specimens from soft tissue such as breast, kidney, liver, lung, thyroid, lymph nodes and various soft tissue masses performed under image guidance techniques for soft tissue biopsy.

CONTRAINDICATIONS

None known

WARNINGS (W)

⚠ **W1.** Good medical judgement should be exercised in considering biopsy on patients who are receiving anticoagulant therapy or who have bleeding disorders. Use in these patient populations may result in severe bleeding.

⚠ **W2.** The instrument is intended to provide breast tissue for histologic examination with partial or complete removal of the imaged abnormality. The extent of histologic abnormality cannot be reliably determined from its mammographic appearance. Therefore, the extent of removal of the imaged evidence of an abnormality does not predict the extent of removal of a histologic abnormality e.g., malignancy. When the sampled abnormality is not histologically benign, it is essential that the tissue margins be examined for completeness of removal using standard surgical procedures.

⚠ **W3.** If the needle are damaged or bent, do not use. Attempt to correct bent parts may cause breakage. Do not use excessive force against resistance, it may damage or fracture the AprioCore plus and result in materials remaining in the body.

⚠ **W4.** AprioCore plus is intended for use in image-guided procedures only. Verify the device location using imaging techniques after every device adjustment. Do not use AprioCore plus in Magnetic Resonance Imaging (MRI) equipment. The device is magnetic and exposure to MRI may cause unintended device movement.

⚠ **W5.** One potential complication of coaxial image-guided biopsy and unintentional activation of biopsy device is inadvertent perforation of adjacent structures. Inadvertent perforation may result in site specific complications which may consist of bleeding, including bleeding from the large artery, hematoma, hemorrhage, infection, pain, cardiac tamponade, pneumothorax, hemothorax and hemothorax.

PRECAUTIONS (P)

⚠ **P1.** The device must be used by a physician who is completely familiar with the indications, contraindications, limitations, typical findings and possible side effects of core needle biopsy, in particular, those relating to the specific organ being biopsied.

⚠ **P2.** Before use together with other biopsy instruments, check the compatibility in length and gauge together with AprioCore plus (see table 1).

⚠ **P3.** Do not test fire AprioCore plus into the air before a procedure. Test fire without support of tissue may cause damage to the needle/cannula tip.

⚠ **P4.** After use, the components may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

DIRECTION FOR USE

Preparation:

1. If a coaxial needle is used, determine which version is preferred. Suitable coaxial needle is described in table 1.
2. Use aseptic technique while removing AprioCore plus from the package.
3. Before use, inspect the device. If the needle is damaged or bent, do not use. ⚠ **W1, W2, W3**

Biopsy Procedure:

1. If a coaxial needle is used, advance the coaxial needle to the area to be biopsied as per the coaxial needle Instruction for Use. Confirm that the needle is in position using image guidance. ⚠ **W4**

2. Prior to insertion of AprioCore plus, load the device by pulling the trigger button until a firm click is felt, indicating that the device is ready for use (figure 1).
3. Confirm that the stylet is fully retracted and that the specimen release button is in the rearmost position (figure 2). The specimen notch shall be completely covered by the cutting cannula.
4. Advance the needle to the area to be biopsied (figure 3). Do not advance the slotted stylet until AprioCore plus is in position. **⚠ W3, W4, W5**
5. While maintaining needle position, advance the slotted stylet with thumb to expose the specimen notch within the area to be biopsied (figure 4). **⚠ W3, W4, W5**
6. Confirm that the specimen release button is in its retracted position (figure 5). Press the trigger button past the first stop to capture the tissue specimen.
7. Withdraw the needle.
8. To remove the tissue specimen, push the specimen release button forward to expose the specimen notch (figure 6).
9. For multiple specimens, repeat step 5-11. **⚠ W3**

OTHER INFORMATION

Serious incident that has occurred in relation to the device shall be reported in accordance with national laws and regulations. Serious incidents shall be reported to AprioMed AB (legal manufacturer) and to the appropriate regulatory agency such as the competent authority of the Member State (EU/EEA) in which the user/patient is established.

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AprioCore® is a registered trademark of AprioMed AB. Patent Pending
Mosaic Surgical Ltd