Corbo Coaxial[®] Needle

Coaxial biopsy needle


- (EN) Instruction for use
- (ES) Instrucciones de uso
- (SV) Bruksanvisning
- (IT) Istruzioni per l'uso
- (NL) Gebruiksaanwijzing
- (DE) Gebrauchsanweisung
- (FR) Mode d'emploi
- (NO) Bruksanvisning

AprioMed

(EN) Table 1. (ES) Tabla 1. (SV) Tabell 1. (IT) Tabell 1. (NL) Tabel 1. (DE) Tabelle 1. (FR) Tableau 1. (NO) Tabell 1.

(EN) REF. (ES) REF. (SV) REF. (IT) RIF. (NL) REF. (DE) REF. (FR) RÉF. (NO) REF.	(EN) Gauge x Total length (ES) Calibre x Longitud total (SV) Gauge x Total längd (IT) Calibro x Lunghezza totale (NL) Dikte x Totale lengte (DE) Dicke x Gesamtlänge (FR) Jauge x Longueur totale (NO) Diameter x Total lengde	(EN) Free length (ES) Longitud libre (SV) Fri längd (IT) Lunghezza libera (NL) Vrije lengte (DE) Freie Länge (FR) Longueur libre (NO) Fri lengde	(EN) Coaxial cannula ID (ES) D. i. de la cánula coaxial (SV) Koaxialkanyl ID (IT) Diametro int. cannula coassiale (NL) ID coaxiale canule (DE) ID der koaxialen Kanüle (FR) DI de la canule coaxiale (NO) Koaksial kanyle ID	(EN) Suitable biopsy device size (Gauge/Length) (ES) Tamaño adecuado del dispositivo de biopsia (Calibre/Longitud) (SV) Passande storlek på biopsienhet (Gauge/Längd) (IT) Dimensioni idonee del dispositivo per biopsia (Calibro/Lunghezza) (NL) Geschikte grootte van het biopsie-instrument (Dikte/Lengte) (DE) Geeignete Größe des Biopsiegeräts (Dicke/Länge) (FR) Taille appropriée de l'instrument de biopsie (Jauge/Longueur) (NO) Størrelse på egnet biopsienhet (Diameter/Lengde)
 CCN17G-7	17G x 7 cm	4.6 cm	1.33 mm	18G ≥ 8.5 cm
 CCN17G-12	17G x 12 cm	9.6 cm	1.33 mm	18G ≥ 13.5 cm
 CCN17G-17	17G x 17 cm	14.6 cm	1.33 mm	18G ≥ 18.5 cm
 CCN19G-7	19G x 7 cm	4.6 cm	0.93 mm	20G ≥ 8.5 cm
 CCN19G-12	19G x 12 cm	9.6 cm	0.93 mm	20G ≥ 13.5 cm
 CCN19G-17	19G x 17 cm	14.6 cm	0.93 mm	20G ≥ 18.5 cm

(EN)

1. Corbo Coaxial® Needle
2. End Cap
3. 22G Anesthesia needle
4. Hollow trocar stylet
5. Coaxial cannula
6. Total length
7. Free length

(SV)

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

ENGLISH

DEVICE DESCRIPTION

Corbo Coaxial® Needle is a manually operated, sterile, single use, coaxial needle with an integrated 22G anesthesia needle. Corbo Coaxial® Needle can be used with biopsy instruments of different sizes, see Table 1. Description of Corbo Coaxial® Needle parts, see Table 2.

INDICATION FOR USE

As a guiding needle in obtaining core biopsy samples from soft tissue such as liver, kidney, spleen, lymph nodes and various soft tissue lesions.

CONTRAINDICATIONS

Not intended for use in bone.

WARNINGS (W)

- ⚠ W1. Good medical judgement should be exercised in considering biopsy on patients who are receiving anticoagulant therapy or who have bleeding disorders. Use in these patient populations may result in severe bleeding.
- ⚠ W2. Do not use excessive force against resistance. If the needle is damaged or bent, do not use. Retract the bent part simultaneously with the rest of the device to prevent needle breakage and materials remaining in the body.
- ⚠ W3. Corbo Coaxial® Needle is intended for use in image-guided procedures only. Verify the device location using imaging techniques after every device adjustment. Do not use Corbo Coaxial® Needle in Magnetic Resonance Imaging (MRI) equipment. The device is magnetic and exposure to MRI may cause unintended device movement.
- ⚠ W4. One potential complication of coaxial image-guided biopsy is inadvertent perforation of adjacent structures. Inadvertent perforation may result in site specific complications which may consist of bleeding, including bleeding from the large artery, hematoma, hemorrhage, infection, pain, cardiac tamponade, pneumothorax, hemothorax, and hemothorax.

PRECAUTIONS (P)

- ⚠ P1. Corbo Coaxial® Needle must be used by a physician who is completely familiar with the indications, contraindications, limitations, typical findings and possible side effects of core needle biopsy, in particular, those relating to the specific organ being biopsied.
- ⚠ P2. Before use together with other equipment, check the compatibility in length and gauge together with Corbo Coaxial® Needle.
- ⚠ P3. Before use, inspect the device for damaged needle point, bent shaft or other imperfections that would prevent proper function. If the needle components are damaged or bent or do not move together easily, do not use.
- ⚠ P4. Administration of local anesthesia must be conducted by, or under supervision of physician with experience of regional anesthesia in accordance with the drug and syringe manufacturer instructions.
- ⚠ P5. After use, the components may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

DIRECTION FOR USE

Preparation

The biopsy procedure must be performed using appropriate aseptic technique.

Before use together with other biopsy instruments check the compatibility in length and gauge together with Corbo Coaxial® Needle, see Table 1. Confirm that the distal end of the biopsy instrument extends through the coaxial cannula.

Before use, inspect the device for damaged needle point, bent shaft or other imperfections that would prevent proper function. If the needle components are damaged, bent or do not move together easily, do not use.

Biopsy Procedure (⚠ W1-W4, ⚠ P1-P5)

1. Remove end cap, attach syringe with anesthesia to the 22G anesthesia needle Luer. Insert Corbo Coaxial® Needle superficially and inject local anesthetics at entry site (**Figure 1**).
2. Advance Corbo Coaxial® Needle towards the target anatomy using image guided technique (**Figure 2**). NOTE! The 22G anesthesia needle protrudes 30mm out of the coaxial cannula. Administer anesthesia through the 22G anesthesia needle during advancement. Administration of local anesthesia must be conducted in accordance with drug and syringe manufacturer instructions.
3. When the target anatomy is reached, confirm correct position using image guidance technique.
4. Disconnect the 22G anesthesia needle, advance the hollow trocar stylet together with the coaxial cannula over the 22G anesthesia needle to the target lesion so that the etch-marking on the 22G anesthesia needle is visible (**Figure 3**).
5. Disconnect the hollow trocar stylet and remove it together with the 22G anesthesia needle (**Figure 4**).
6. Perform biopsy through the coaxial access (**Figure 5**).
7. If repositioning of the coaxial cannula is needed, reinsert the hollow trocar stylet to support the coaxial cannula during repositioning. Confirm correct position using image guidance technique.
8. When the procedure is completed, reinsert the hollow trocar stylet and remove the coaxial cannula (**Figure 6**).

Other information

Serious incident that has occurred in relation to the device shall be reported in accordance with national laws and regulations. Serious incidents shall be reported to AprioMed AB (legal manufacturer) and to the appropriate regulatory agency such as the competent authority of the Member State (EU/EEA) in which the user/patient is established.

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Rx Only

(EN) Caution: Federal (U.S.A) law restricts this device to sale by or on the order of a physician.

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(EN) Single use. Do not re-use. After use, the product performance and characteristics may degrade, and can cause infection in patient.

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(EN) Do not re-sterilize. Re-sterilization may have impact on product characteristics and performance.

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(EN) Do not use if package is damaged. Device is delivered sterile. Product in damaged package may be contaminated and can cause infection in patient.

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(EN) MR unsafe. Do not use device in Magnetic Resonance Imaging (MRI) equipment. The device is magnetic and exposure to MRI may cause unintended device movement.

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Corbo Coaxial® Needle is a registered trademark of AprioMed AB.
Patent Pending.

UK Responsible Person:

Mosaic Surgical Ltd
9 York Way
High Wycombe
HP12 3PY
United Kingdom

STERILE EO

Rx Only

**UK
CA
0086**

AprioMed AB
Virdings Allé 28
SE-754 50 Uppsala, Sweden
Tel: +46 (0) 18 430 14 40
Fax: +46 (0) 770 33 99 72
info@apriomed.com
www.apriomed.com

AprioMed