
SeeStar®

Guiding Device

- (EN) Instruction for use
- (SV) Bruksanvisning
- (NL) Handleiding
- (FR) Mode d'emploi
- (ES) Instrucciones de uso
- (IT) Istruzioni per l'uso
- (DE) Gebrauchsanweisung
- (NO) Bruksanvisning

AprioMed

(EN) Figures (SV) Figurer (NL) Afbeeldingen (FR) Figure (ES) Figuras (IT) Figura (DE) Abbildungen (NO) Figurer


Fig. 1

Fig. 2


Fig. 3

(EN) Table 1. (SV) Tabell 1. (NL) Tabel 1. (FR) Tableau 1. (ES) Tabla 1. (IT) Tabella 1. (DE) Tabelle 1. (NO) Tabell 1.

REF	Gauge	Color	Inner diameter	Recommended needle size
13312	12G	Blue 	2,9 mm 0.114 inch	12G 13G
13314	14G	Green 	2,2 mm 0.087 inch	14G 15G
13317	17G	Pink 	1,6 mm 0.063 inch	17G 18G
13320	20G	Yellow 	1,0 mm 0.039 inch	20G 21G

ENGLISH

INDICATION FOR USE

A passive guide for guided invasive access to tissue that is to be examined and/or treated.

WARNINGS (W)

- ⚠ **W1.** If using SeeStar® in conjunction with other equipment, first check that the equipment is compatible with SeeStar®. See table for maximum permissible needle diameter.
- ⚠ **W2.** Note that the needle holder is 4 cm long, which must be taken into account when calculating the puncture length of the needle (Fig.3).

PRECAUTIONS (P)

- ⚠ **P1.** SeeStar® should only be used by, or under the supervision of a physician.
- ⚠ **P2.** Once SeeStar® has been taped to the skin, it could be dislodged if the skin moves. If the patient has been moved, a control image should be taken to check its location.
- ⚠ **P3.** If an external needle is used, the needle holder should be loosened after puncture so as not to restrict unnecessarily the body's normal internal movements
- ⚠ **P4.** After use, the components may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

DIRECTIONS FOR USE

SeeStar® can be used to simplify and improve the precision of the regular procedure for CT guided puncture.

- Apply local anaesthetic.
- Make a small skin incision prior to the puncture.
- The entry point is determined and marked along the laser line. SeeStar is applied with the tip of the needle guide directly over the entry point and the hub open (Fig.1)
- Make sure the needle guide is in the image plane by centering/aligning the proximal end of the metal needle guide with the laser line (Fig.1). The linear artefact visible in the CT image marks the puncture route through the body/ the intended needle path. (Fig.2)
- The angle of the needle guide is adjusted in order to reach the target and the needle guide is locked in the desired position by tightening the hub.
- If repeat punctures are required from the same area, the SeeStar® can be left in its position. Other areas can also be punctured from the same entry point by angling the needle guide.

OTHER INFORMATION

Serious incident that has occurred in relation to the device shall be reported in accordance with national laws and regulations. Serious incidents shall be reported to AprioMed AB (legal manufacturer) and to the appropriate regulatory agency such as the competent authority of the Member State (EU/EEA) in which the user/patient is established.

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Rx Only

(EN) Caution: Federal (U.S.A) law restricts this device to sale by or on the order of a physician.

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(EN) Single use. Do not re-use. After use, the product performance and characteristics may degrade, and can cause infection in patient.

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(EN) Do not re-sterilize. Re-sterilization may have impact on product characteristics and performance.

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(EN) Do not use if package is damaged. Device is delivered sterile. Product in damaged package may be contaminated and can cause infection in patient.

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(EN) MR unsafe. Do not use AprioCore plus in Magnetic Resonance Imaging (MRI) equipment. The device is magnetic and exposure to MRI may cause unintended device movement.

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SeeStar® is a registered trademark of AprioMed AB.
Mosaic Surgical Ltd