


Bonopty[®]

Bone Biopsy System, 12 gauge



Optimize tactility and efficiency during
your procedures with Bonopty[®] products

AprioMed

Bonopty® Bone Biopsy System, 12 gauge

- New, larger gauge gains access into the bone, even through thick cortical bone
- Coaxial system facilitates multiple sampling, or treatment of lesions
- Achieve excellent core samples with few crush artifacts
- Successfully sample somewhat sclerotic lesions
- Secure sample in biopsy cannula during retrieval using Core Lock



Use of Bonopty® Penetration Set and Bonopty® Biopsy Set to sample a lytic lesion with bone structures in lumbar vertebrae.



Before using Bonopty® read the Instructions for Use which accompany the product. Bonopty® is a registered trademark of AprioMed AB, Uppsala, Sweden. Bonopty® is manufactured under one or more of the following patents: US 5,423,824; US 5,810,826; US 5,868,684; EP 0 794 733; JP 3710475; SE 470 177. Other patents pending.

ORDERING INFORMATION

Bonopty® Bone Biopsy System, 12G

Product number	Product Description	Minimum Order Quantity
12-1272	Bonopty® Penetration Set, 12 gauge, 10 cm	5 units
12-1273	Bonopty® Biopsy Set, 13 gauge, 16.5 cm	5 units
12-1274	Bonopty® Extended Drill, 13 gauge, 16.5 cm	5 units

Bonopty® Bone Biopsy System, 14G

Product number	Product Description	Minimum Order Quantity
10-1062	Bonopty® Penetration Set, 14 gauge, 6.5 cm	5 units
10-1072	Bonopty® Penetration Set, 14 gauge, 9.5 cm	5 units
10-1082	Bonopty® Penetration Set, 14 gauge, 13 cm	5 units
10-1063	Bonopty® Biopsy Set, 15 gauge, 13 cm	5 units
10-1073	Bonopty® Biopsy Set, 15 gauge, 16 cm	5 units
10-1083	Bonopty® Biopsy Set, 15 gauge, 19.5 cm	5 units
10-1074	Bonopty® Extended Drill, 15 gauge, 16 cm	5 units