

DuoGuide® II

GUIDEWIRE

TaperTip™

DILATOR

Optimize tactility and efficiency during the first stage of the nephrostomy procedure



Indication for use

DuoGuide® II Guidewire: For use in angiographic procedures to introduce and position catheters and interventional devices within the coronary and peripheral vasculature.

TaperTip™ Dilator: For use in percutaneous nephrostomy and transhepatic cholangiography procedures.



AprioMed

DUOGUIDE®II GUIDEWIRE AND TAPERTIP™ DILATOR


Designed to make the first stage of the nephrostomy procedure simpler and less painful

DuoGuide®II Guidewire

- Catheter support without compromising formability
- Achieve accurate placement through visualization
- Introduce through a 0.9 mm needle


TaperTip™ Dilator

- Achieve dilation in a single step
- Lessen pain with the long tapered tip
- Reduce risk of penetrating the collecting system



ORDERING INFORMATION

Product number	Product description	Minimum Order Quantity
13002	DuoGuide®II Guidewire Tip Flexibility Std super-elastic with shapeable tip <ul style="list-style-type: none">• Ø 0.46 mm / 0.018"• 900 mm length• 45 mm soft platinum tip• Nitinol shaft	5 units
13003	DuoGuide®II Guidewire Tip Flexibility High super-elastic with shapeable tip <ul style="list-style-type: none">• Ø 0.46 mm / 0.018"• 900 mm length• 80 mm soft platinum tip• Nitinol shaft	5 units
13100	TaperTip™ Dilator radioopaque, compatible with 0.018" guidewire <ul style="list-style-type: none">• Inner Ø 0.51 mm• Outer Ø 3.3 mm• 30 mm tapered tip• 240 mm length	10 units

 Prior to using DuoGuide®II or TaperTip™, read the Instructions for Use accompanying the products. DuoGuide®II and TaperTip™ are not available for sale in the United States. DuoGuide®II and TaperTip™ are trademarks of AprioMed AB, Uppsala, Sweden.

HEAD OFFICE

AprioMed AB
Virdings Allé 28
SE-754 50 Uppsala, Sweden
Phone: +46 18 430 14 40
Free Fax (Int'l): +46 770 339972

US SALES OFFICE

AprioMed, Inc.
2 Palmer Dr, Suite 1
Londonderry, NH 03053, USA
Phone: +1 603 421 0875
Fax (US toll free): 866 239 1797



info@apriomed.com / www.apriomed.com

APIRIOMED develops, manufactures, markets and sells innovative medical devices and related services within the field of interventional radiology. We aim to deliver, in close collaboration with healthcare practitioners, innovative quality tools to achieve the most favourable solutions for radiologists worldwide. AprioMed is headquartered in Uppsala, Sweden, with a sales subsidiary in Londonderry, NH, USA. Our products are sold in more than 40 countries by AprioMed (direct sales in US, UK, Ireland and Sweden) and our local distributors.